

## Pearl White **Professional Teeth Whitening System**

Order #: UDPWPRO



### **Pearl White Stand Alone Features:**

**Easy to Use Hand Held Devices**  
**Blue LED Enhances the Performance of the Whitening Gel**  
**Advanced Design Allows Constant Light Without Generating Heat Providing Fast performance**  
**Optimum Treatment Time 30 to 40 Minutes**  
**Option of Using 1 or 2 Mouth Pieces Simultaneously**  
**Four Treatment Timing Selections**  
**Medical Grade Power Supply**  
**Accepts Worldwide AC Current Automatically**  
**Alleviate Teeth Sensitivity**  
**No Eye Protection Required**  
**Power Mains – 100-240 VAC, 50/60 Hz**  
**Made in U.S.A.**

### **ACCESSORIES INCLUDED WITH THE PURCHASE OF THE EQUIPMENT**

- One Power Cord
- One Teeth Whitening LED Mouth Piece , (1 years Warranty)
- 50 Gel Syringes
- 25 Waiting Room Brochures
- One Manual & Warranty Card (3 years Warranty), *Accessories Warranty 90 Days*

### **Benefits:**

Provide your client beauty inside and out. The Carbamide Peroxide gel is based on a tested and true formulation made of high quality ingredients. The gel rapidly produces a noticeably brighter and whiter smile. Most people will experience 4-6 shades lighter.

Teeth whitening can be a positive and less time consuming experience for you and your client. A small amount of whitening gel is placed on the mouth piece in the LED location. The client positions the mouth piece in their mouth up against their teeth and the technician activates the light simple and effective.



Before



After